

## Freedom® 600 bag

Your technician will supply you with a purpose-made bag to carry your Freedom® 600 cylinder and regulator.

1. Make sure the cylinder is secured inside the bag using the Velcro straps attached (see picture below).
2. Push the trumpet end of your cannula through the access hole in the bag and attach the cannula to the regulator.
3. Once the cylinder has been secured in the bag, fasten the zip fully.
4. You can monitor the amount of oxygen in the cylinder while it is inside the bag through the window provided.
5. Use one of the three handles attached to carry the bag.



### Caution

When using the Freedom® 600, use the bag to carry the cylinder at all times.

When transporting in a pushchair or car, make sure the cylinder bag is safely secured at all times.