

## Fig.1 The 3 main types of concentrators

On average concentrators weigh around 25kg (55.1lb) and are 50–70cm (19.5 – 27.5ins) high. They make a slight humming noise, like a fridge, when switched on.

### AirSep Newlife Elite concentrator

- A On/Off switch
- B Flow meter/Adjustment knob
- C Oxygen outlet
- D Indicator lights
- E Handle grip
- F Filter
- G Power cord
- H Circuit breaker reset button
- I Hour meter



### DeVilbiss concentrator

- A On/Off switch
- B Flow meter/Adjustment knob
- C Oxygen outlet
- D Indicator lights
- E Handle grip
- F Filter
- G Power socket
- H Circuit breaker reset button



### Respironics Millennium concentrator

- A On/Off switch
- B Flow meter/Adjustment knob
- C Oxygen outlet
- D Indicator lights
- E Handle grip
- F Filter
- G Power socket
- H Hour meter

